



NOVAIR



ANALYTICAL GRADE NITROGEN AND COMPRESSED AIR SOLUTIONS FOR LIFE SCIENCES

814 WURLITZER DR., NORTH TONAWANDA, NY 14120 USA

716-564-5165

INFO@NOVAIR-USA.COM



WWW.NOVAIR-USA.COM



Known For the Best Service in the Industry

GREAT USER EXPERIENCE – WE AT NOVAIR ARE COMMITTED TO MAKING SURE WE ARE ALWAYS THERE FOR YOU, WHENEVER AND WHEREVER YOU NEED OUR NOVAIR SKILLED TECHNICIANS, TO MAKE SURE YOUR FACILITY’S EQUIPMENT IS ALWAYS RUNNING SMOOTHLY WITH MINIMAL DOWNTIME.

LOCAL SERVICE TECHS -- NOVAIR’S OWN NATIONWIDE TEAM OF FACTORY-TRAINED FIELD SERVICE ENGINEER SPECIALISTS ARE THE BEST IN THE BIZ AND COME STRAIGHT FROM THE SOURCE


GREAT PREVENTATIVE MAINTENANCE OPTIONS – THE NOVAIR-CARE PREVENTATIVE SERVICE PLAN OPTIONS INCLUDE 3 DIFFERENT LEVELS TO CHOOSE FROM, OR A ‘CHOOSE YOUR OWN’ PLAN AS WELL, SO THAT YOU RECEIVE REGULAR MAINTENANCE IN THE MANNER YOU DESIRE, TO KEEP YOUR EQUIPMENT RUNNING SMOOTH FOR YEARS AND YEARS TO COME.



WWW.NOVAIR-USA.COM



FOR OVER 45 YEARS, THE NOVAIR BRAND HAS BEEN RENOWNED GLOBALLY FOR MANUFACTURING RELIABLE AND INNOVATIVE MEDICAL SYSTEMS FOR ALL LIFE SCIENCE AND HOSPITAL GAS DISTRIBUTION AND GENERATION NEEDS. AT NOVAIR USA, WE TAKE PRIDE IN MANUFACTURING STATE-OF-THE-ART OXYGEN, NITROGEN, MEDICAL AIR AND VACUUM GENERATOR SYSTEMS RIGHT HERE IN THE UNITED STATES. COMMITTED TO EXCELLENCE, OUR SYSTEMS ARE METICULOUSLY DESIGNED AND MANUFACTURED IN STRICT ACCORDANCE WITH NFPA99 STANDARDS TO ENSURE THE HIGHEST LEVEL OF SAFETY AND RELIABILITY.



NOVAIR's superior Nitrogen Generators and turnkey Nitrogen Generation Systems with tailored solutions for your life sciences facility translates to the highest quality laboratory gas on demand

FEATURES:

- **LEADING TECHNOLOGY** – STATE-OF-THE-ART MODULARITY, ALLOWING ANYONE IN YOUR FACILITY TO EASILY UPGRADE THE UNIT AND EXPAND YOUR NITROGEN SUPPLY AS YOUR FACILITY EXPANDS
- **HIGH-QUALITY** – GET THE HIGHEST NITROGEN CONCENTRATION ALONG WITH OIL-FREE, CLEAN, DRY COMPRESSED AIR IN ALL NITROGEN GENERATING SYSTEMS
- **ENERGY EFFICIENT** – REDUCE OPERATIONAL COSTS BY ONLY GETTING THE NITROGEN YOU NEED, WHEN YOU NEED IT, WITH EASY-TO-USE TOUCHSCREENS AND CONSTANT MONITORING FOR CONSISTENT FLOW WITH NO WASTE
- **EQUISITE DESIGN** – REGARDED AS THE INDUSTRY LEADING IN-HOUSE PSA NITROGEN GENERATOR, THE VERSATILE, DURABLE, AND EXPERTLY ENGINEERED NITROSWING® PROVIDES MAXIMUM VALUE WITH A MINIMUM FOOTPRINT
- **QUIET** – WHISPER QUIET OPERATING TECH MEANS NEAR NOISE-FREE

BENEFITS:

- **HAVE GAS ALWAYS, ON-DEMAND, AND REST EASY** KNOWING YOU HAVE CONTINUOUS FLOW OF THE HIGHEST-QUALITY NITROGEN THAT YOUR FACILITY NEEDS, BUT AT A FRACTION OF THE COST OF OUTSIDE SUPPLY
- **TREMENDOUS ROI** BY GENERATING YOUR OWN NITROGEN FROM AMBIENT AIR FOR THE PRICE OF ONLY ELECTRICITY TO RUN THE GENERATOR
- **SAVE BUDGET TO SPEND ON OTHER NEEDS** BY ENDING RELIANCE ON COSTLY DELIVERY WHILE ALSO REDUCING YOUR CARBON FOOTPRINT
- **GAIN PEACE-OF-MIND** AS THE NITROSWING® CONSTANTLY OPERATES AT PEAK PERFORMANCE WITH REDUCED MAINTENANCE AND LITTLE TO NO DOWNTIME, SO YOU CAN FOCUS ON WHAT MATTERS MOST
- **FEEL CONFIDENT** IN YOUR PURCHASE OF QUALITY USA AND EUROPEAN DESIGNED AND MANUFACTURED EQUIPMENT BY THE INDUSTRY LEADER WITH A PROVEN TRACK RECORD



NOVAIR's NITROSWING patented true modular PSA Nitrogen Generator turnkey solution

WITH THE NITROSWING, NOT ONLY DO YOU GET THE HIGHEST QUALITY NITROGEN GENERATOR ON THE MARKET THAT FILLS TONS OF CYLINDERS FOR THE PRICE OF ELECTRICITY, EASILY PLUGGING RIGHT INTO YOUR FACILITY'S PIPING SYSTEM, BUT YOU ALSO GET THE SYSTEM YOU WANT, THE WAY YOU WANTED, GROWING AS YOU GROW. A MORE RELIABLE SOLUTION COMPARED WITH MEMBRANE TECHNOLOGY, NITROSWING'S TIME-TESTED PSA TECHNOLOGY, TOGETHER WITH NOVAIR'S RIGOROUS PRODUCT TESTING STANDARDS, MEANS YOU CAN BREATHE EASY EASY KNOWING THE MACHINE WILL CONSISTENLY OPERATE AT PEAK PERFORMANCE ALL. ALL NITROSWING GENERATORS COME STANDARD WITH A SYSTEM THAT IS TRULY MODULAR, ALLOWING FOR A RAPID INCREASE IN OXYGEN PRODUCTION SHOULD YOUR EXPANDING LABORATORY OR CLINIC NEED IT. UNLIKE THE CLAIMS OF OUR COMPETITORS, THE NITROSWING CAN BE VERY EASILY UPGRADED BY ADDING EXTRA COMPONENTS WITH VIRTUALLY NO ADJUSTMENTS NEEDED, MEANING ANY PERSON ON YOUR TEAM CAN DO IT.

- THE HIGHEST QUALITY MODULAR PSA NITROGEN GENERATOR ON THE MARKET, NOVAIR'S NITROSWING BOASTS SUPERIOR TECHNOLOGY, SLEEK ENGINEERING DESIGN AND A SMALL FOOTPRINT, ALL BACKED BY THE BEST SERVICE OUT THERE
- EASILY INTEGRATED WITH YOUR FACILITY'S INFRASTRUCTURE
- NOVAIR'S NITROSWING'S PSA TECHNOLOGY PROVIDES THE HIGHEST PURITY NITROGEN AVAILABLE
- THE NITROSWING IS THE ONLY TRULY MODULAR EASILY EXPANDABLE N2 PSA GENERATING SYSTEM ON THE MARKET -- UNLIKE THE CLAIMS OF OUR COMPETITORS, THE NITROSWING CAN BE EASILY UPGRADED BY ANY TEAM MEMBER IN YOUR FACILITY



Why choose NOVAIR - USA?

- **MADE IN THE USA AND EUROPE:** OUR SYSTEMS ARE PROUDLY DESIGNED AND MANUFACTURED IN ITALY, FRANCE, AND THE UNITED STATES, ENSURING THE HIGHEST GLOBAL QUALITY WITH ADHERENCE TO LOCAL STANDARDS.
- **APPROVALS:** MEDICAL AIR AND VACUUM SYSTEMS MANUFACTURED IN OUR NOVAIR USA PRODUCTION FACILITIES ARE DESIGNED AND MANUFACTURED IN STRICT ACCORDANCE WITH NFPA99 GUIDELINES, ENSURING STRICT REGULATORY COMPLIANCE.
- **EXCEPTIONAL SUPPORT:** WE OFFER COMPREHENSIVE CUSTOMER SUPPORT, INCLUDING INSTALLATION, MAINTENANCE, AND TECHNICAL ASSISTANCE, TO KEEP YOUR FACILITY RUNNING SMOOTHLY
- **CUSTOM SOLUTIONS:** NOVAIR-USA CAN TAILOR OUR SYSTEMS TO MEET YOUR SPECIFIC REQUIREMENTS, ENSURING OPTIMAL PERFORMANCE AND EFFICIENCY
- **CHOOSE NOVAIR-USA FOR TOP-OF-THE-LINE OXYGEN, NITROGEN, MEDICAL AIR AND VACUUM SYSTEMS THAT PRIORITIZE SAFETY, EFFICIENCY, AND RELIABILITY;** CONTACT US TODAY TO DISCUSS HOW OUR CUTTING-EDGE TECHNOLOGIES CAN BENEFIT YOUR HEALTHCARE OR LABORATORY FACILITY



WWW.NOVAIR-USA.COM

